

We are '*Professionals in Process/ Product Technologies*' with following expertise-

- **Conceptualization and Engineering of Projects**
- **Module and skid building**
- **Custom-built solutions**
- **Innovative hygiene Products for Pharma/ Dairy Industry**
- **Expertise in Pharma, Food, Natural resource based projects**
- **Business partners of world class products from Europe and USA**

Our Offerings-

- **Bio-Reactors And Fermenters**
- **Bio - Waste Decontamination Systems / Sterilization Systems / Inactivation Systems**
- **PWDS, WFIDS, CSDS Systems**
- **CIP Plants**
- **Pharmaceutical Machinery and Bio-pharma Custom Built Modules**
- **Mixing And Blending Vessels**
- **Dairy/Starch Projects**
- **Drying And Evaporation Projects**
- **Biorefineries**
- **Piping Engineering**
- **Plug-in Temperature / Brix / pH / Pressure/ Flow Control Modules**
- **Electrical And Automation Engineering**
- **Flow Control Equipment**
- **All Types Of Valves And Valve Automation**

Application

VISILAMP is high quality Sight & Light glass assembly in concealed SS304 enclosure that enables process personnel to have visual feel of the process in following types of equipment

- Bioreactors, Fermenters
- Mixing Vessels
- Blending Vessels
- Storage Tanks, Silos
- Chemical reaction vessels
- Chemical equipment
- Calendrias / Condensers
- Vapour Separators



Features

- Can be installed as sight / light glass
- Low heat output compared with light energy
- Low Power requirement, <15 watt
- Designed for hygienic applications
- Luminous intensity for spot visibility
- Spring loaded gasket
- Easy accessibility for the lamp assembly
- Concealed SS 304 enclosure, can open only by bottom lock nut

Technical Specifications

MOC of Enclosure	SS 304
External Finish	Satin <0.8 Ra
Lamp Type	13W PLC Havells
MOC of Gasket	Silicone Rubber
Glass	Borosilicate, 5mm thick, toughened.

FOR MORE DETAILS CONTACT

Office & Works :-SM BIOSYSTEMS
S. No. 28 / 27, Dhayari, Narhe - Dhayari Road,
Behind Oneness Controls, Dhayari,
Pune - 411041.
Maharashtra INDIA

Tel.-020 64701545, 020 24690018
Fax-020 24690655
Email- mktg@smbiosystems.com
response@smbiosystems.com
Visit us at www.smbiosystems.com