

We are '*Professionals in Process/ Product Technologies*' with following expertise-

- **Conceptualization and Engineering of Projects**
- **Module and skid building**
- **Custom-built solutions**
- **Innovative hygiene Products for Pharma/ Dairy Industry**
- **Expertise in Pharma, Food, Natural resource based projects**
- **Business partners of world class products from Europe and USA**

**Our Offerings-**

- **Bio-Reactors And Fermenters**
- **Bio - Waste Decontamination Systems / Sterilization Systems / Inactivation Systems**
- **PWDS, WFIDS, CSDS Systems**
- **CIP Plants**
- **Pharmaceutical Machinery and Bio-pharma Custom Built Modules**
- **Mixing And Blending Vessels**
- **Dairy/Starch Projects**
- **Drying And Evaporation Projects**
- **Biorefineries**
- **Piping Engineering**
- **Plug-in Temperature / Brix / pH / Pressure/ Flow Control Modules**
- **Electrical And Automation Engineering**
- **Flow Control Equipment**
- **All Types Of Valves And Valve Automation**

## PURIFLUX- DW/PW/WFI Distribution System



### Description & Applications-

PURIFLUX enables consumption of Water-for-injection/ Purified/ Distilled water with continuous flow of water at the consumption point as per cGMP norms. PURIFLUX fitting is specially designed to be welded in the main loop. This is to avoid running main loop up to the point of use and to avoid using bigger valves. The unwanted dead legs are avoided. This system is terminated at consumption point using specially fabricated 3 way SS316L body diaphragm valve.



### Features-

1. Designed as per cGMP norms, Hygienic
2. MOC SS316L, Electropolished
3. Main loop not required to run up to consumption point
4. Various sizes available as per consumption
5. No welded parts
6. Cost efficient

FOR MORE DETAILS CONTACT

Office & Works :-SM BIOSYSTEMS  
S. No. 28 / 27, Dhayari, Narhe - Dhayari Road,  
Behind Oneness Controls, Dhayari,  
Pune - 411041.  
Maharashtra INDIA

Tel.-020 64701545, 020 24690018  
Fax-020 24690655  
Email- [mktg@smbiosystems.com](mailto:mktg@smbiosystems.com)  
[response@smbiosystems.com](mailto:response@smbiosystems.com)  
Visit us at [www.smbiosystems.com](http://www.smbiosystems.com)